

# Senior Health Monitoring & Caring



- Health monitoring wearables can be used to track medical condition, record data, and provides trends
- An emergency condition may occur without the senior realising it (e.g. when in sleep)






- Consensus can then be obtained from relative to allow a home visit by the medical personnel
- Cost effective solution without the need to set-up monitoring centre with agents to take/make calls








## Layer 1 - Product Infrastructure




- Health monitoring wearable in the form of a watch or chest strap
- 2. Storage of health data collected from the wearable
- Computation to generate reports and trends
- 4. Phone registry of medical doctors
- 5. Phone book of relatives

#### Layer 2 - Sensors




- GPS that provides location data in case the senior citizen is not at the registered address
- 2. Temperature sensor to detect body temperature
- Heat beat sensor to detect the cardiac condition of the senior citizen
- There could be other sensors involved that provides other health data

## Layer 3 - Connectivity



- The wearable connects to the cloud using Bluetooth connection to a smart phone
- 2. The smart phone can then relay the data to the cloud through Wi-Fi
- 3. The doctor will be informed by text message through a mobile phone network
- 4. The relative will be informed by text message through a mobile phone network

## Level 4 - Analytics



- 1. Graphs are produced for the collected data
- 2. Highlights of abnormal conditions against the senior citizen's historical data and that of the community will be provided
- Any detected trending in the health data will be highlighted

#### Level 5 – Smart Apps



- Smart phone app to view health data records and graphs
- 2. Smart phone app to inform the doctor and the relatives about abnormal medical condition
- 3. Consensus by the relative to pay the doctor for a home visit
- 4. Report of home visit to be sent to relatives